

Histone H2B Mouse Monoclonal Antibody(Mix-mA)

Catalog No.: RA10012

Basic Information

Information

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	14KD

Applications

Recommended Dilution


WB 1:1,000-3,000

IHC 1:200-500


Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)Hela, 2)3T3 , 3)Raw264.7, 4)Rat Brain, 5)Rat Kidney with Histone H2B(Mix-mATM) Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human breast carcinoma using Histone H2B Mouse mAb diluted at 1:500.

Background

Histone H2B is one of the 5 main histone proteins involved in the structure of chromatin in eukaryotic cells.