

# Cystatin C Mouse Monoclonal Antibody(5F1)

Catalog No.: RA10117

## Basic Information

### Information

Reactivity	H
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	14KD

### Applications


### Recommended Dilution

WB	1:1,000-2,000
IHC	1:100-200
ELISA	1:10,000-20,000


## Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Immunohistochemical analysis of paraffin-embedded Human Kidney Tissue using Cystatin C Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Brain Tissue using Cystatin C Mouse mAb diluted at 1:200.