

Cystatin C Mouse Monoclonal Antibody(3B12)

Catalog No.: RA10118

Basic Information

Information

Reactivity	H
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	14KD

Applications

Recommended Dilution

WB	1:1,000-2,000
ELISA	1:10,000-20,000

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice