

# Luciferase Mouse Monoclonal Antibody(2G11)

Catalog No.: RA10142

## Basic Information

### Information

Reactivity	Firefly
Immunogen	Synthetic Peptide
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	60KD

### Applications


### Recommended Dilution

WB	1:1,000-2,000
----	---------------

### Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Western blot analysis of Firefly Luciferase Transfected Hela Cell Lysate using Luciferase Mouse mAb diluted at 1:2000.