

# Cytokeratin 6 Mouse Monoclonal Antibody(4A1)

Catalog No.: RA10358

## Basic Information

### Information

Reactivity	H
Immunogen	Synthetic Peptide
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	56-60KD

### Applications


### Recommended Dilution

WB	1:1,000
IHC	1:100-200

## Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Western blot analysis of 1)HeLa, 2)HepG2 ,with Cytokeratin Mouse mAb diluted at 1:2,000.

## Background

Cytokeratins are proteins of keratin-containing intermediate filaments found in the intracytoplasmic cytoskeleton of epithelial tissue.