

EFHD1 Mouse Monoclonal Antibody(3G2)

Catalog No.: RA10363

Basic Information

Information

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	27KD

Applications


Recommended Dilution

WB	1:2,000
IF	1:100-200

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1) Mouse spleen tissue, 2) Rat spleen tissue, IF analysis of HeLa with EFHD1 Mouse mAb (Left) with EFHD1 Mouse mAb diluted at 1:3,000 and DAPI (Right) diluted at 1:100.

Background

EFHD1 is an EF-hand domain-containing protein that displays increased expression during neuronal differentiation.