

Histone H3(mono methyl K79) Mouse Monoclonal Antibody(1C2)

Catalog No.: RA10376

Basic Information

Information

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	15KD

Applications


Recommended Dilution

WB	1:1,000-3,000
IHC	1:200-500


Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)Rat Testis tissue, 2)Raw264.7, 3)Mouse Brain tissue with Histone H3 (mono methyl K79) Mouse MAb diluted at 1:2000.



Immunohistochemical analysis of paraffin-embedded Human Skin using Histone H3(mono methyl K79) Mouse mAb diluted at 1:500.

Background

Histone H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.