

TurboGFP Mouse mAb(Mix-mA)**Catalog No.: RA10410****Basic Information****Information**

Reactivity	Species independent
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	~26KD


Applications**Recommended Dilution**

WB	1:20,000-100,000
-----------	------------------

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of Recombinant MBP-TurboGFP Protein using TurboGFP Mouse Monoclonal Antibody diluted at 1:50,000.

Background

TurboGFP is an improved variant of the green fluorescent protein CopGFP cloned from copepod *Pontellina plumata*.