

**Rabbit IgG(Fc specific,Heavy Chain) Mouse Monoclonal
Antibody(1E1)
Catalog No.: RA10428**

Basic Information

Information

Reactivity	Rabbit
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	55KD

Applications


Recommended Dilution

WB	1:5,000-20,000
IHC	1:200-500

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of Reduced Rabbit IgG, H+L with Rabbit IgG(Fc specific, Heavy Chain) Mouse Monoclonal Antibody diluted at 1:10,000.