

Histone H3(Mono Methyl Arg2) Rabbit Polyclonal Antibody

Catalog No.: RA20089

Basic Information

Information

Reactivity	M
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	15-17KD

Applications


Recommended Dilution

WB	1:500-1,000
----	-------------

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)Raw264.7 with Histone H3 (Mono Methyl Arg2)Rabbit pAb diluted at 1:1,000.

Background

Histone H3 is one of the five main histone proteins involved in the structure of chromatin in eukaryotic cells. Core component of nucleosome. Nucleosomes wrap and compact DNA into chromatin, limiting DNA accessibility to the cellular machineries which require DNA as a template. Histones thereby play a central role in transcription regulation, DNA repair, DNA replication and chromosomal stability.