

GSK 3 β Rabbit Polyclonal Antibody(F41、 F43)

Catalog No.: RA20126

Basic Information

Information

Reactivity	M,R
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	46KD

Applications


Recommended Dilution

WB	1:1,000-2,000
IHC	1:200-500


Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)Mouse Brain Tissue, 2)Mouse Heart Tissue, 3)Rat Brain Tissue, 4)Rat Heart Tissue with GSK 3 β Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Mouse Lung Tissue using GSK 3 β Rabbit pAb diluted at 1:500.

Background

Glycogen synthase kinase 3 beta, also known as GSK3B, is an enzyme that in humans is encoded by the GSK3B gene.