

## EPG5 Rabbit Polyclonal Antibody(121-240)

Catalog No.: RA20128

### Basic Information

#### Information

Reactivity	H,M,R
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	170-292KD

#### Applications


#### Recommended Dilution

WB	1:1,000-2,000
IHC	1:200-500


#### Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Western blot analysis of Rat Liver Tissue with EPG5 Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Mouse Liver Tissue using EPG5 Rabbit pAb diluted at 1:500.

## Background

N/A.