

CALHM1 Rabbit Polyclonal Antibody(A191)

Catalog No.: RA20138

Basic Information

Information

Reactivity	H,R
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	40KD

Applications



Recommended Dilution

WB	1:500-1,000
IHC	1:200-500

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1) Human Brain Tissue, 2) Rat Brain Tissue with CALHM1 Rabbit pAb diluted at 1:2,000.

Immunohistochemical analysis of paraffin-embedded Human Hepatocarcinoma using CALHM1 Rabbit pAb diluted at 1:500.

Background

Calcium homeostasis modulator 1 (CALHM1) is a pore-forming membrane protein that in humans is encoded by the CALHM1 gene.