

MEK2 Rabbit Polyclonal Antibody(F189)

Catalog No.: RA20166

Basic Information

Information

Reactivity	H,M,R
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	45KD

Applications



Recommended Dilution

WB	1:1,000-2,000
IHC	1:200-500

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)HeLa, 2)3T3, 3)PC12 with MEK2 Rabbit pAb diluted at 1:2,000.

Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma using MEK2 Rabbit pAb diluted at 1:500.

Background

MEK1 and MEK2, also called MAPK or Erk kinases, are dual-specificity protein kinases that function in a mitogen activated protein kinase cascade controlling cell growth and differentiation. Catalyzes the concomitant phosphorylation of a threonine and a tyrosine residue in a Thr-Glu-Tyr sequence located in MAP kinases. Activates ERK1 and ERK2 MAP kinases.