

# AMPK $\beta$ 2 Rabbit Polyclonal Antibody(F110)

Catalog No.: RA20177

## Basic Information

### Information

Reactivity	H,M,R
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	30KD

### Applications


### Recommended Dilution

WB	1:1,000-2,000
----	---------------

## Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Western blot analysis of 1)293T, 2)HepG2, 3)Mouse Heart Tissue, 4)Rat Heart Tissue with AMPK $\beta$ 2 Rabbit pAb diluted at 1:2,000.

## Background

5' AMP-activated protein kinase or AMPK or 5' adenosine monophosphate-activated protein kinase is an enzyme that plays a role in cellular energy homeostasis.