

# STAT5b Rabbit Polyclonal Antibody(F255)

Catalog No.: RA20193

## Basic Information

### Information

Reactivity	H,M,R
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	90KD

### Applications


### Recommended Dilution

WB	1:1,000-2,000
IHC	1:100-200


## Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Western blot analysis of 1)K562, 2)Hela, 3)3T3, 4)Mouse Heart Tissue, 5)PC12, 6)Rat Heart Tissue with STAT5bRabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Spleen Tissue using STAT5b Rabbit pAb diluted at 1:200.

## Background

Signal transducer and activator of transcription 5B is a protein that in humans is encoded by the STAT5B gene. STAT5B orthologs have been identified in most placentals for which complete genome data are available.