

**c-Jun Rabbit Polyclonal Antibody(F127)****Catalog No.: RA20194****Basic Information****Information**

<b>Reactivity</b>	H,M
<b>Immunogen</b>	Recombinant Protein
<b>Host</b>	Rabbit
<b>Isotype</b>	IgG
<b>Storage Buffer &amp; Condition</b>	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
<b>Observed MW</b>	43KD

**Applications****Recommended Dilution**



<b>WB</b>	1:1,000-2,000
<b>IHC</b>	1:50-100

**Preparation & Storage**

<b>Storage</b>	Store at -20°C. Stable for one year from the date of shipment.
<b>Shipping</b>	Bule Ice

## Experimental Data

错误！未指定文件名。



Western blot analysis of 1)293T, 2)Hela, 3)3T3 with c-Jun Rabbit pAb diluted at 1:2,000.

Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using c-Jun Rabbit pAb diluted at 1:200.

## Background

c-Jun is a member of the Jun Family containing c-Jun, JunB and JunD, and is a component of the transcription factor AP-1 (activator protein-1).