


ASIC1 Rabbit Polyclonal Antibody(A180)**Catalog No.: RA20203****Basic Information****Information**

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	80KD


Applications**Recommended Dilution****WB** 1:1,000-2,000**IHC** 1:50-100**Preparation & Storage**

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)293T, 2)Mouse Brain Tissue, 3)Rat Brain Tissue with ASIC1 Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using ASIC1 Rabbit pAb diluted at 1:200.

Background

Amiloride-sensitive cation channel 2,neuronal is a protein that in humans is encoded by the ACCN2 gene.