

CCKBR Rabbit Polyclonal Antibody(A214)

Catalog No.: RA20219

Basic Information

Information

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	50KD

Applications


Recommended Dilution

WB	1:1,000-2,000
----	---------------

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)Mouse Brain Tissue, 2)Rat Brain Tissue, 3)Human Brain Tissue, with CCKBR Rabbit pAb diluted at 1:2,000.

Background

Receptor for gastrin and cholecystokinin. The CCK-B receptors occur throughout the central nervous system where they modulate anxiety,analgesia,arousal,and neuroleptic activity.