

KCNN3(SK3) Rabbit Polyclonal Antibody(A246)

Catalog No.: RA20244

Basic Information

Information

Reactivity	H,R
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	82KD

Applications


Recommended Dilution

IHC	1:100-200
-----	-----------

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using KCNN3(SK3)Rabbit pAb diluted at 1:200.

Background

Forms a voltage-independent potassium channel activated by intracellular calcium. Activation is followed by membrane hyperpolarization.