

# KV $\beta$ 1 Rabbit Polyclonal Antibody(A255)

Catalog No.: RA20250

## Basic Information

### Information

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	50KD

### Applications


### Recommended Dilution

IHC	1:100-200
-----	-----------

## Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Kvb1 Rabbit pAb diluted at 1:200.

## Background

N/A