

## KV1.5 Rabbit Polyclonal Antibody(A257)

Catalog No.: RA20252

### Basic Information

#### Information

Reactivity	H
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	67KD

#### Applications


#### Recommended Dilution

WB	1:1,000-2,000
----	---------------

#### Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Western blot analysis of 293T with KV1.5 Rabbit pAb diluted at 1:2,000.