

KV1.8 Rabbit Polyclonal Antibody(A258)

Catalog No.: RA20253

Basic Information

Information

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	58KD

Applications


Recommended Dilution

WB	1:1,000-2,000
IHC	1:100-200


Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)Rat Brain Tissue, 2)Mouse Brain Tissue with Kv1.8 Rabbit pAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using Kv1.8 Rabbit pAb diluted at 1:200.

Background

N/A