

KV11.3 Rabbit Polyclonal Antibody(A263)

Catalog No.: RA20257

Basic Information

Information

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	83135KD

Applications



Recommended Dilution

WB	1:1,000-2,000
IHC	1:100-200

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)Rat Brain Tissue, 2)Mouse Brain Tissue , 3)HeLa with KV11.3 Rabbit pAb diluted at 1:2,000

Immunohistochemical analysis of paraffin-embedded Rat BrainTissue using Kv11.3 Rabbit pAb diluted at 1:200.

Background

N/A