

DUSP6 Rabbit Polyclonal Antibody

Catalog No.: RA20258

Basic Information

Information

Reactivity	H,M,R
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	42,44KD

Applications


Recommended Dilution

WB	1:1,000-5,000
----	---------------

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)HepG2, 2)3T3 , 3)Rat Heart Tissue with DUSP6 Rabbit pAb diluted at 1:3,000.