

# JAK1 Rabbit Polyclonal Antibody(F172)

Catalog No.: RA20286

## Basic Information

### Information

Reactivity	H,M,R
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	130KD

### Applications


### Recommended Dilution

IHC	1:100-200
-----	-----------

## Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using JAK1 Rabbit pAb diluted at 1:200.

## Background

Members of the Janus family of tyrosine kinases (Jak1, Jak2, Jak3 and Tyk2) are activated by ligands binding to a number of associated cytokine receptors.