

Branchiostoma belcheri PAX9 Rabbit Polyclonal Antibody

Catalog No.: RA20330

Basic Information

Information

Reactivity	Branchiostoma blecheri (H,M,R)
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	38KD

Applications


Recommended Dilution

WB	1:500-1,000
IF	1:100-200

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 1)PC-3 Cell, 2)HepG2 Cell, 3) Jurkat Cell Lysate using PAX9 Rabbit pAb diluted at 1:1,000.

Background

Transcription factor required for normal development of thymus, parathyroid glands, ultimobranchial bodies, teeth, skeletal elements of skull and larynx as well as distal limbs.