

# Syndecan-4 Rabbit Polyclonal Antibody

## Catalog No.: RA20344

### Basic Information

#### Information

Reactivity	H,M,R
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	26KD

#### Applications


#### Recommended Dilution

WB	1:1,000-2,000
----	---------------

#### Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



Western blot analysis of Mouse Liver Tissue Lysate using Syndecan-4 Rabbit Polyclonal Antibody diluted at 1:2,000.

## Background

The protein encoded by this gene is a transmembrane (type I) heparan sulfate proteoglycan that functions as a receptor in intracellular signaling. The encoded protein is found as a homodimer and is a member of the syndecan proteoglycan family.