

GAPDH Rabbit Polyclonal Antibody

Catalog No.: RCA59

Basic Information

Information

Reactivity	H,M,R,Mk,Ch,Hm,Sh,Xenopus
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	37KD

Applications


Recommended Dilution

WB	1:5,000
IHC	1:200

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of 293T (1), Rat brain (2), NIH 3T3(3), Sheep Muscle(4),Rabbit testis(5) with GAPDH rabbit polyclonal diluted at 1:20,000.

Background

Glyceraldehyde 3 phosphate dehydrogenase (GAPDH) is well known as one of the key enzymes involved in glycolysis. GAPDH is constitutively expressed in almost all tissues at high levels, therefore antibodies against GAPDH are useful as loading controls for Western Blotting. Some physiological factors, such as hypoxia and diabetes, increase GAPDH expression in certain cell types.