

SRT-Tag Mouse Monoclonal Antibody(4H9)

Catalog No.: RTA19

Basic Information

Information

Reactivity	N/A
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	N/A

Applications


Recommended Dilution

WB	1:5,000
----	---------

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of SRT-Recombinant protein with diluted at 1)1:5,000 2)1:10,000

Background

The SRT-tag is a peptide epitope comprised of 10 amino acids(TFIGAIATDT) derived from the crystalline surface layer protein of *Rickettsia typhi*.