

RFP-Tag Mouse Monoclonal Antibody(Mix-mA™)

Catalog No.: RTA20

Basic Information

Information

Reactivity	N/A
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	N/A

Applications

Recommended Dilution

WB	1:5,000
IHC	1:200
IP	1:200

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of recombinant RFP protein with diluted at 1:10,000

Background

The red fluorescent protein cloned from *Discosoma* coral (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the *Aequorea* jellyfish.