

S1-Tag Mouse Monoclonal Antibody(7B10)

Catalog No.: RTA37

Basic Information

Information

Reactivity	ALL
Immunogen	Synthetic Peptide
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	N/A

Applications

Recommended Dilution

WB	1:5,000-10,000
----	----------------

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of Recombinant S1-Tag Protein with S1-Tag Mouse mAb diluted at 1:5,000.

Background

S1 tag is an epitope tag composed of a nine-residue peptide, NANNPDWDF, derived from the hepatitis B virus preS1 region.