

CBP Tag Mouse Monoclonal Antibody(Mix-mA™)

Catalog No.: RTA42

Basic Information

Information

Reactivity	N/A
Immunogen	Recombinant Protein
Host	Mouse
Isotype	IgG1
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	N/A

Applications


Recommended Dilution

WB	1:5,000
IP	1:200

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of CBP recombinant protein with diluted at 1) 1:5,000 2) 1:10,000

Background

Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.