

# RFP-Tag Rabbit Polyclonal Antibody

## Catalog No.: RTA58

### Basic Information

#### Information

Reactivity	N/A
Immunogen	Recombinant Protein
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	N/A

#### Applications


#### Recommended Dilution

WB	1:5,000
----	---------

#### Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

## Experimental Data



2ug RFP fusion protein+ Primary antibody dilution at 1、 1:5,000 2、 1:10,000

## Background

The red fluorescent protein cloned from *Discosoma coral* (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the *Aequorea jellyfish*.