

SBP-tag Rabbit Polyclonal Antibody

Catalog No.: RTA60

Basic Information

Information

Reactivity	ALL
Immunogen	Synthetic Peptide
Host	Rabbit
Isotype	IgG
Storage Buffer & Condition	1mg/ml in PBS, pH 7.4, containing 0.02% sodium azide and 50% glycerol.
Observed MW	ALL

Applications

Recommended Dilution

WB	1:5,000-10,000
----	----------------

Preparation & Storage

Storage	Store at -20°C. Stable for one year from the date of shipment.
Shipping	Bule Ice

Experimental Data



Western blot analysis of Recombinant SBP-Tag protein using SBP-Tag Rabbit pAb diluted at 1:20,000
Western blot analysis of Recombinant SBP-Tag protein using SBP-Tag Rabbit pAb diluted at 1:20,000

Background

Streptavidin Binding Peptide (SBP) tag is very versatile and can be inserted into expression vectors for easy detection of recombinant proteins and purification by use of affinity columns.